

CLYBL Antibody



PACO08476

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CLYBL

Uniprot

Q8N0X4

Synonyms:

citrate lyase beta like; CLYBL; CLB; bA134O15.1 ;

Immunogen:

Human CLYBL.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A