

CNOT10 Antibody



PACO08491

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

CNOT10

Uniprot

Q9H9A5

Synonyms:

CCR4-NOT transcription complex, subunit 10; CNOT10; DKFZp434K115; FLJ12890; FLJ13165 ;

Immunogen:

Human CNOT10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A