

COL25A1 Antibody



PACO08510

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

COL25A1

Uniprot

Q9BXS0

Synonyms:

51kd or 64kd;

Immunogen:

Human COL25A1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A