

# COLEC11 Antibody



PACO08517

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat, Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COLEC11

**Uniprot**

Q9BWP8

**Synonyms:**

collectin sub-family member 11; COLEC11; CL-K1-I; CL-K1-II; CL-K1-IIa; CL-K1-IIb; CLK1; DKFZp686N1868; MGC3279 ;

**Immunogen:**

Human COLEC11.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A