

COPS3 Antibody



PACO08525

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

COPS3

Uniprot

Q9UNS2

Synonyms:

COP9 constitutive photomorphogenic homolog subunit 3 (Arabidopsis); COPS3; CSN3; SGN3 ;

Immunogen:

Human COPS3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A