

COQ10A Antibody



PACO08531

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

COQ10A

Uniprot

Q96MF6

Synonyms:

coenzyme Q10 homolog A (S. cerevisiae); COQ10A; FLJ32452 ;

Immunogen:

Human COQ10A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A