

COX15 Antibody



PACO08542

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

COX15

Uniprot

Q7KZN9

Synonyms:

COX15 homolog, cytochrome c oxidase assembly protein (yeast);

Immunogen:

Human COX15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A