

# COX18 Antibody



PACO08545

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COX18

**Uniprot**

Q8N8Q8

**Synonyms:**

COX18 cytochrome c oxidase assembly homolog (S. cerevisiae); COX18; COX18HS; FLJ38991; MGC126733 ;

**Immunogen:**

Human COX18.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A