

COX5A Antibody



PACO08547

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

COX5A

Uniprot

P20674

Synonyms:

cytochrome c oxidase subunit Va; COX5A; COX; COX-VA; VA ;

Immunogen:

Human COX5A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A