

# COX6B2 Antibody



PACO08551

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

COX6B2

**Uniprot**

Q6YFQ2

**Synonyms:**

cytochrome c oxidase subunit VIb polypeptide 2 (testis); COX6B2; COXVIB2; CT59; MGC119094 ;

**Immunogen:**

Human COX6B2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A