

CPT1A Antibody



PACO08581

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

CPT1A

Uniprot

P50416

Synonyms:

carnitine palmitoyltransferase 1A (liver); CPT1A; CPT1; CPT1-L; L-CPT1 ;

Immunogen:

Human CPT1A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A