

CRABP2 Antibody



PACO08585

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CRABP2

Uniprot

P29373

Synonyms:

cellular retinoic acid, binding protein 2; CRABP2; CRABP-II; RBP6 ;

Immunogen:

Human CRABP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A