

CRADD Antibody



PACO08586

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CRADD

Uniprot

P78560

Synonyms:

CASP2 and RIPK1 domain containing adaptor with death domain; CRADD; MGC9163; RAIDD ;

Immunogen:

Human CRADD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A