

CREB1 Antibody



PACO08590

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CREB1

Uniprot

P16220

Synonyms:

cAMP responsive element binding protein 1; CREB1; CREB; MGC9284 ;

Immunogen:

Human CREB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A