

CREBL2 Antibody



PACO08595

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CREBL2

Uniprot

O60519

Synonyms:

cAMP responsive element binding protein-like 2; CREBL2; MGC117311; MGC138362 ;

Immunogen:

Human CREBL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A