

CREBZF Antibody



PACO08596

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CREBZF

Uniprot

Q9NS37

Synonyms:

CREB/ATF bZIP transcription factor; CREBZF; SMILE; ZF ; ; Zhangfei

Immunogen:

Human CREBZF.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A