

PACO08601

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CRIP1

Uniprot

Q9P021

Synonyms:

cysteine-rich PDZ-binding protein; CRIP1; HSPC139 ;

Immunogen:

Human CRIP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

Product Images

N/A

N/A