

CRISP3 Antibody



PACO08603

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CRISP3

Uniprot

P54108

Synonyms:

cysteine-rich secretory protein 3; CRISP3; Aeg2; CRISP-3; CRS3; MGC126588; SGP28; dJ442L6.3 ;

Immunogen:

Human CRISP3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A