

CRK Antibody



PACO08604

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CRK

Uniprot

P46108

Synonyms:

v-crk sarcoma virus CT10 oncogene homolog (avian); CRK; CRKII ;

Immunogen:

Human CRK.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A