

CRLS1 Antibody



PACO08605

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CRLS1

Uniprot

Q9UJA2

Synonyms:

cardiolipin synthase 1; CRLS1; C20orf155; CLS1; GCD10; dj967N21.6 ;

Immunogen:

Human CRLS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A