

CRTAP Antibody



PACO08610

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CRTAP

Uniprot

O75718

Synonyms:

cartilage associated protein; CRTAP; CASP; LEPREL3 ;

Immunogen:

Human CRTAP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A