

CSF2RA Antibody



PACO08625

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CSF2RA

Uniprot

P15509

Synonyms:

colony stimulating factor 2 receptor, alpha, low-affinity (granulocyte-macrophage); CSF2RA; CD116; CDw116; CSF2R; CSF2RAX; CSF2RAY; CSF2RX; CSF2RY; GM-CSF-R-alpha; GMCSFR; GMR; MGC3848; MGC4838 ;

Immunogen:

Human CSF2RA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A