

# CSGALNACT2 Antibody



PACO08626

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CSGALNACT2

**Uniprot**

Q8N6G5

**Synonyms:**

chondroitin sulfate N-acetylgalactosaminyltransferase 2; CSGALNACT2; CSGaINAcT-2; DKFZp686H13226; FLJ43310; GALNACT-2; GALNACT2; MGC40204; PRO0082 ;

**Immunogen:**

Human CSGALNACT2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A