

CSH2 Antibody



PACO08628

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CSH2

Uniprot

P0DML3

Synonyms:

chorionic somatomammotropin hormone 2; CSH2; CS-2; CSB; hCS-B ;

Immunogen:

Human CSH2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A