

# CSNK2A1 Antibody



PACO08635

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CSNK2A1

**Uniprot**

P68400

**Synonyms:**

casein kinase 2, alpha 1 polypeptide; CSNK2A1; CK2A1; CKII ;

**Immunogen:**

Human CSNK2A1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A