

# CST9L Antibody



PACO08645

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CST9L

**Uniprot**

Q9H4G1

**Synonyms:**

cystatin 9-like; CST9L; FLJ92394; bA218C14.1 ;

**Immunogen:**

Human CST9L.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A