

CTDSP2 Antibody



PACO08655

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CTDSP2

Uniprot

O14595

Synonyms:

CTD (carboxy-terminal domain, RNA polymerase II, polypeptide A) small phosphatase 2; CTDSP2; OS4; PSR2; SCP2 ;

Immunogen:

Human CTDSP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A