

# CTNNA1 Antibody



PACO08658

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

CTNNA1

**Uniprot**

P35221

**Synonyms:**

catenin (cadherin-associated protein), alpha 1, 102kDa; CTNNA1; CAP102; FLJ36832; FLJ52416 ;

**Immunogen:**

Human CTNNA1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A