

CUX1 Antibody



PACO08681

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CUX1

Uniprot

Q13948

Synonyms:

cut-like homeobox 1; CUX1; CASP; CDP; CDP/Cut; CDP1; COY1; CUTL1; CUX; Clox; Cux/CDP; GOLIM6; Nbla10317; p100; p110; p200; p75 ;

Immunogen:

Human CUX1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A