

CXCL16 Antibody



PACO08689

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CXCL16

Uniprot

Q9H2A7

Synonyms:

chemokine (C-X-C motif) ligand 16; CXCL16; CXCLG16; SR-PSOX; SRPSOX ;

Immunogen:

Human CXCL16.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A