

CXCR4 Antibody



PACO08692

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CXCR4

Uniprot

P61073

Synonyms:

chemokine (C-X-C motif) receptor 4; CXCR4; CD184; D2S201E; FB22; HM89; HSY3RR; LAP3; LCR1; LESTR; NPY3R; NPYR; NPYRL; NPYY3R; WHIM ; FUSIN; CXCR4B

Immunogen:

Human CXCR4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A