

CYB5B Antibody



PACO08696

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CYB5B

Uniprot

O43169

Synonyms:

cytochrome b5 type B (outer mitochondrial membrane); CYB5B; CYB5-M; CYPB5M; DKFZp686M0619; OMB5 ;

Immunogen:

Human CYB5B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A