

CYB5R4 Antibody



PACO08699

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CYB5R4

Uniprot

Q7L1T6

Synonyms:

cytochrome b5 reductase 4; CYB5R4; NCB5OR; RP4-676J13.1; dJ676J13.1 ;

Immunogen:

Human CYB5R4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A