

CYP2C8 Antibody



PACO08714

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CYP2C8

Uniprot

P10632

Synonyms:

cytochrome P450, family 2, subfamily C, polypeptide 8; CYP2C8; CPC8; MP-12/MP-20 ;

Immunogen:

Human CYP2C8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A