

CYP2J2 Antibody



PACO08716

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CYP2J2

Uniprot

P51589

Synonyms:

cytochrome P450, family 2, subfamily J, polypeptide 2; CYP2J2; CPJ2 ;

Immunogen:

Human CYP2J2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A