

CYP46A1 Antibody



PACO08721

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CYP46A1

Uniprot

Q9Y6A2

Synonyms:

cytochrome P450, family 46, subfamily A, polypeptide 1; CYP46A1; CP46; CYP46 ;

Immunogen:

Human CYP46A1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A