

# CYTL1 Antibody



PACO08730

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CYTL1

**Uniprot**

Q9NRR1

**Synonyms:**

cytokine-like 1; CYTL1; C17; C4orf4 ;

**Immunogen:**

Human CYTL1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A