

DAAM1 Antibody



PACO08732

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

DAAM1

Uniprot

Q9Y4D1

Synonyms:

dishevelled associated activator of morphogenesis 1; DAAM1; FLJ41657; KIAA0666 ;

Immunogen:

Human DAAM1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A