

DAB2 Antibody



PACO08733

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DAB2

Uniprot

P98082

Synonyms:

disabled homolog 2, mitogen-responsive phosphoprotein (Drosophila); DAB2; DOC-2; DOC2; FLJ26626 ;

Immunogen:

Human DAB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A