

DAP3 Antibody



PACO08739

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DAP3

Uniprot

P51398

Synonyms:

death associated protein 3; DAP3; DAP-3; DKFZp686G12159; MGC126058; MGC126059; MRP-S29; MRPS29; bMRP-10 ;

Immunogen:

Human DAP3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A