

DARS2 Antibody



PACO08743

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DARS2

Uniprot

Q6PI48

Synonyms:

aspartyl-tRNA synthetase 2, mitochondrial; DARS2; ASPRS; FLJ10514; LBSL; MT-ASPRS; RP3-383J4.2 ;

Immunogen:

Human DARS2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A