

DAZL Antibody



PACO08746

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DAZL

Uniprot

Q92904

Synonyms:

deleted in azoospermia-like; DAZL; DAZH; DAZL1; DAZLA; MGC26406; SPGYLA ;

Immunogen:

Human DAZL

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A