

DBH Antibody



PACO08749

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DBH

Uniprot

P09172

Synonyms:

dopamine beta-hydroxylase (dopamine beta-monoxygenase); DBH; DBM ;

Immunogen:

Human DBH.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A