

DCBLD2 Antibody



PACO08755

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DCBLD2

Uniprot

Q96PD2

Synonyms:

discoidin, CUB and LCCL domain containing 2; DCBLD2; CLCP1; ESDN ;

Immunogen:

Human DCBLD2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A