

DCTD Antibody



PACO08763

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

DCTD

Uniprot

P32321

Synonyms:

dCMP deaminase; DCTD; MGC111062 ;

Immunogen:

Human DCTD.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A