

DDA1 Antibody



PACO08772

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DDA1

Uniprot

Q9BW61

Synonyms:

DET1 and DDB1 associated 1; DDA1; C19orf58; MGC2594; PCIA1 ;

Immunogen:

Human DDA1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A