

# DDOST Antibody



PACO08780

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DDOST

**Uniprot**

P39656

**Synonyms:**

dolichyl-diphosphooligosaccharide-protein glycosyltransferase; DDOST; AGE-R1; KIAA0115; MGC2191; OK/SW-cl.45; OST; OST48; WBP1 ;

**Immunogen:**

Human DDOST.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A