

DDX20 Antibody



PACO08784

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DDX20

Uniprot

Q9UHI6

Synonyms:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 20; DDX20; DKFZp434H052; DP103; GEMIN3 ;

Immunogen:

Human DDX20.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A