

DDX23 Antibody



PACO08786

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

DDX23

Uniprot

Q9BUQ8

Synonyms:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 23; DDX23; MGC8416; PRPF28; U5-100K; U5-100KD; prp28 ;

Immunogen:

Human DDX23.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A