

# DDX25 Antibody



PACO08788

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DDX25

**Uniprot**

Q9UHL0

**Synonyms:**

DEAD (Asp-Glu-Ala-Asp) box polypeptide 25; DDX25; GRTH ;

**Immunogen:**

Human DDX25.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A